



www.deroyal.com



Content

Product Listing	Page
SURGICAL ACCESSORIES	
Vessel Loops	3
Disposable Light Handle Systems	4
Light Handle Adapters	5
Suture Boots	6
Needle Counter Systems	7-8
BLADES & SCALPELS	
Disposable Safety Scalpels with Retractable Blades	9-10
Swann-Morton® Scalpels, Blades and Handles	11-14
Scalpel Blade Remover	8
DeBlade™ Needle Counters	8
SURGICAL CONVENIENCE	
Instrument Pads	15
Skin Markers and Rulers	16
Surgical Utility Markers	16
Dual Tip Marker	16
WOUND CLOSURE	
Skin Staplers	17-19
Skin Stapler Extractor	19
SURGICAL ESSENTIALS	
Sponge Counter Rack	20
Instrument Pouch	21
Cautery Tip Cleaner	21
Instrument Wipes	21
SAFETY PRODUCTS	
Scalpel Holder	22
Safety Strainers	22
Needle Cap Caddie™	23
HANDS FREE PASSING	
SharpStop™ Sharps Transfer Tray	24
SURGICAL CLIPS AND CLAMPS	
Multi-Purpose Tube Holder	25
Towel Clamp	26
Tube Positioner	26
Breathing Circuit Support Clamp	27
Anesthesia Drape "Peach" Clip	27
DRAPES & COVERS	
C-Arm Drapes	28
Camera Drapes	29
Cassette Drapes	30-31
Band Bags	32
Isolation Drape	33
Foot-Switch Cover	33
Slush Drape	34
Case Cart Cover	34

* All products in this catalog are sterile unless otherwise noted.



The art of surgery requires action and decisiveness, starting with selection of the supplies and instruments. DeRoyal provides surgical accessories that not only get the job done, but let you do it the way you want, using thoughtful designs that offer more convenience at less expense.

Vessel Loops

- High quality medical-grade silicone
- Strong and durable, yet gentle
- Non-porous to slide smoothly and easily over vessels
- Radiopaque element ensures x-ray detectability
- Color-coded to assist in vessel detection
- Dispenser boxes minimize storage
- Three sizes available
- Two 16" strands included in each pack:
 - Minis are 1.3mm wide and 0.9mm thick
 - Maxis are 2.5mm wide and 1.0mm thick
 - Super Maxis are 4.7mm wide and 1.3mm thick
- Neatly wrapped on a plastic card for tangle-free presentation and placement in sterile field
- Not made with natural rubber latex



Product No.	Description	Quantity
30-711	Vessel Loops, White Mini, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-712	Vessel Loops, Red Mini, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-713	Vessel Loops, Blue Mini, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-714	Vessel Loops, Yellow Mini, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-721	Vessel Loops, White Maxi, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-722	Vessel Loops, Red Maxi, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-723	Vessel Loops, Blue Maxi, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-724	Vessel Loops, Yellow Maxi, 2/Card	2/Pkg, 25 Pkg/Bx, 2 Bx/Cs
30-731	Vessel Loops, White Mini, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-732	Vessel Loops, Red Mini, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-733	Vessel Loops, Blue Mini, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-734	Vessel Loops, Yellow Mini, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-741	Vessel Loops, White Maxi, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-742	Vessel Loops, Red Maxi, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-743	Vessel Loops, Blue Maxi, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-744	Vessel Loops, Yellow Maxi, 2/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs
30-1728	Vessel Loops, Blue Super Maxi, 1/Card	2/Pkg, 10 Pkg/Bx, 1 Bx/Cs

FLEXIBLE LIGHT HANDLE SYSTEMS



26-012

SurgiClick™ Flexible Light Handle Covers

- Reduces the cost of processing, sterilizing and storing light handles
- Flat packaging requires less than half the space and waste bulk than other products
- Adapters are provided to fit all brands of surgical lights, promoting standardization
- Patented design snaps firmly and easily into place to minimize contamination risks
- Unique ergonomic design minimizes risk of splitting or tearing of plastic while ensuring a no-slip grip
- Rigid flange allows quick, easy application

Product No.	Description	Quantity
26-011	SurgiClick™ Light Handle Covers, Flexible (Snap-on)	1/Pk, 60 Pks/Bx, 2 Bxs/Cs
26-012	SurgiClick™ Light Handle Covers, Flexible (Snap-on)	2/Pk, 30 Pks/Bx, 2 Bxs/Cs
26-013	SurgiClick™ Light Handle Covers, Flexible (Snap-on)	3/Pk, 25 Pks/Bx, 2 Bxs/Cs



AM3611

Saf-T-Glove™ Flexible Light Handle Covers

- Double angle flange for easier application and aiding in keeping cover in place
- Serpentine rib design for reinforcement of flange while applying cover. This reduces the possibility of contamination
- Price is highly competitive within the market of light handle covers
- Application to universal adapters for standardization of application process, providing an easier conversion process
- Friction-fit provides a secure fit that will minimize risk of cover falling off during procedures so as to not contaminate the sterile field

Product No.	Description	Quantity
AM3611	Saf-T-Glove™ Light Handle Covers, Flexible (Friction-Fit)	1/Pk, 50 Pks/Bx, 2 Bxs/Cs
AM3612	Saf-T-Glove™ Light Handle Covers, Flexible (Friction-Fit)	2/Pk, 50 Pks/Bx, 2 Bxs/Cs
AM3613	Saf-T-Glove™ Light Handle Covers, Flexible (Friction-Fit)	3/Pk, 50 Pks/Bx, 2 Bxs/Cs



Rigid Light Handles

- Rigid flange allows quick aseptic application and will bend inward
- Screw-on design is safe and secure
- Semi-rigid plastic is firm yet comfortable to grasp
- Textured handles ensure no-slip grip
- Unique interlocking design eliminates need to package using latex-containing rubber bands or paper bands which add particulates to air
- Special packaging configuration facilitates sterile presentation and minimizes risk of handles falling off sterile field during presentation
- Universal light handle adapters included for almost every model light
- Not made with natural rubber latex

Product No.	Description	Quantity
26-021	Light Handles, Rigid (Screw-on)	1/Pk, 12 Pks/Bx, 4 Bxs/Cs
26-022	Light Handles, Rigid (Screw-on)	2/Pk, 12 Pks/Bx, 2 Bxs/Cs
26-022B	Light Handles, Rigid (Screw-on)	2/Pk, 72 Pks/Cs
26-023	Light Handles, Rigid (Screw-on)	3/Pk, 8 Pks/Bx, 2 Bxs/Cs

LIGHT HANDLE ADAPTERS



Flexible Light Handle Cover Adapters


Light Handle Adapters

DeRoyal's Rigid and Flexible disposable light handle systems includes our universal light handle adapters. These provided adapters fit almost all brands and models of surgical lights, promoting standardization and efficiency.



Rigid Light Handle Adapters

SURGICAL ACCESSORIES



Suture Boots

- Facilitates grasping and tagging without damage to delicate sutures
- Foam cartridge promotes accountability
- Extended adhesive peel-off backing for secure placement and quick removal
- Available in two sizes for a perfect fit:
 - Standard are .062" inside diameter
 - Large are .093" inside diameter
- Radiopaque element ensures x-ray detectability
- Offered in a variety of colors
- Packaged in glassine bag for easy sterile presentation
- Not made with natural rubber latex


Product No.	Description	Quantity
Sterile	Non-Sterile	
30-701*	30-701NS*	Suture Boots, Standard Yellow
30-702*	30-702NS*	Suture Boots, Standard Yellow
30-706*	N/A	Suture Boots, Standard Yellow
30-703**	30-703NS**	Suture Boots, Large Pink
30-704**	30-704NS**	Suture Boots, Large Pink
30-705*	N/A	Suture Boots, Standard Assorted
		10/Pk, 50 Pks/Cs
		6/Pk, 50 Pks/Cs
		10/Pk, 20 Pks/Cs
		10/Pk, 25 Pks/Cs
		6/Pk, 25 Pks/Cs
		10/Pk, 10 Pks/Bx, 2 Bxs/Cs

⚠️ *WARNING - This product can expose you to chemicals including di-isodecyl phthalate, which is known to the State of California to cause birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

⚠️ **WARNING - This product can expose you to chemicals including diisononyl phthalate, which is known to the State of California to cause cancer, and di-isodecyl phthalate, which is known to the State of California to cause birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.



NEEDLE COUNTERS



Deluxe Needle Counters

- Full-surface magnet provides extra space for odd-shaped sharps, pins, blades or needles
- Numbers on magnet allow easy sharps placement and positioning
- Easy-to-read numbers and description facilitate correct count of blades and needles
- See-through cover available, allowing clear box visualization when closed
- Red color signifies biohazard waste
- Not made with natural rubber latex

Product No. Description

Quantity

Magnet

25-0101	Needle Counter, Double Magnet (20 Count 40 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-0104	Needle Counter, Double Magnet (30 Count 60 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-0104W	Needle Counter, Double Magnet (30 Count 60 Capacity) "White Magnet"	50 Pks/Bx, 2 Bxs/Cs
25-0108	Needle Counter, Double Magnet (60 Count 120 Capacity)	27 Pks/Bx, 4 Bxs/Cs
25-0108W	Needle Counter, Double Magnet (60 Count 120 Capacity) "White Magnet"	27 Pks/Bx, 4 Bxs/Cs
25-0201	Needle Counter, Magnet/Adhesive (10 Count 20 Capacity)	50 Pks/Bx, 2 Bxs/Cs

Magnet/Foam Combo

25-0401	Needle Counter, Magnetic/Foam (20 Count 40 Capacity)	20 Pks/Bx, 4 Bxs/Cs
25-0402	Needle Counter, Magnetic/Foam (10 Count 20 Capacity)	25 Pks/Bx, 4 Bxs/Cs
25-0403	Needle Counter, Magnetic/Foam (40 Count 80 Capacity)	12 Pks/Bx, 4 Bxs/Cs
25-0701	Needle Counter, Foam Strips/Magnetic (15 Count 30 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-0702	Needle Counter, Magnetic/Foam (30 Count 60 Capacity)	27 Pks/Bx, 4 Bxs/Cs

Foam

25-0301	Needle Counter, Double Foam (80 Count 160 Capacity)	12 Pks/Bx, 4 Bxs/Cs
25-0501	Needle Counter, Foam/Adhesive (10 Count 20 Capacity)	25 Pks/Bx, 4 Bxs/Cs
25-0502	Needle Counter, Foam/Adhesive (20 Count 40 Capacity)	20 Pks/Bx, 4 Bxs/Cs
25-0504	Needle Counter, Foam/Adhesive (40 Count 80 Capacity)	12 Pks/Bx, 4 Bxs/Cs
25-0801	Needle Counter, Foam Strips (30 Count 60 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-0803	Needle Counter, Foam Strips (60 Count 120 Capacity)	27 Pks/Bx, 4 Bxs/Cs

Clear Top

25-0103	Needle Counter, Single Magnet/Clear Top (10 Count 20 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-0105	Needle Counter, Single Magnet/Clear Top (30 Count 60 Capacity)	27 Pks/Bx, 4 Bxs/Cs
25-0802	Needle Counter, Foam Strips (15 Count 30 Capacity)	50 Pks/Bx, 2 Bxs/Cs

NEEDLE COUNTERS

USING DEROYAL'S DEBLADE REMOVER: THREE EASY STEPS



Step 1: Insert blade into DeBlade mechanism.



Step 2: Slide forward until blade "catches."



Step 3: Slide backward until blade drops safely into box.

DeBlade™ Needle Counters

Product No. Description

Quantity

Magnet

25-1104	DeBlade™ Needle Counter, Double Magnet (30 Count 60 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-1108	DeBlade™ Needle Counter, Double Magnet (60 Count 120 Capacity)	27 Pks/Bx, 4 Bxs/Cs

Foam/Magnet Combo

25-1401	DeBlade™ Needle Counter, Foam Block/Magnetic (20 Count 40 Capacity)	20 Pks/Bx, 4 Bxs/Cs
25-1403	DeBlade™ Needle Counter, Foam Block/Magnetic (40 Count 80 Capacity)	12 Pks/Bx, 4 Bxs/Cs
25-1407	DeBlade™ Needle Counter, Foam Block/Magnetic (70 Count 140 Capacity)	12 Pks/Bx, 4 Bxs/Cs

Foam

25-1502	DeBlade™ Needle Counter, Foam Block/Adhesive (20 Count 40 Capacity)	20 Pks/Bx, 4 Bxs/Cs
25-1504	DeBlade™ Needle Counter, Foam Block/Adhesive (40 Count 80 Capacity)	12 Pks/Bx, 4 Bxs/Cs
25-1701	DeBlade™ Needle Counter, Foam Strips/Magnetic (15 Count 30 Capacity)	50 Pks/Bx, 2 Bxs/Cs
25-1702	DeBlade™ Needle Counter, Foam Strips/Magnetic (30 Count 60 Capacity)	27 Pks/Bx, 4 Bxs/Cs
25-1801	DeBlade™ Needle Counter, Foam Strips (30 Count 60 Capacity)	25 Pks/Bx, 4 Bxs/Cs

Scalpel Blade Remover

- Safely removes blade from metal scalpel handle and contains it for disposal in a sharps container
- Place on the mayo stand or back table in open position for fast, safe blade removal
- Hinged plastic sheath securely locks and contains blade
- Works on any and all scalpel blades and handles
- Single use device
- Not made with natural rubber latex

Product No. Description

Quantity

D6200	Blade Remover, Sterile	50/Cs
D6200NS	Blade Remover, Non-Sterile	100/Cs



BLADES & SCALPELS



Disposable Safety Scalpels With Retractable Blades


- The industry's only patented* safety scalpel that allows a controlled retraction of the entire blade
- Safe and easy to open with temporary forward locking to prevent accidental retractions
- Permanent lock position for disposal now standard on all versions
- 16 different stainless steel blade sizes with color-coded sliders for easy identification
- Ribbed, slip-resistant finger grips for a secure hold in gloved hands, accommodating a variety of hand positions
- High-visibility red color denotes a biohazard to ensure proper disposal in sharps container
- Full blade visibility at all times to reduce the chance of accidental injury
- Designed by Dr. Shackleford, a cardiovascular surgeon, in conjunction with a safety engineer
- Clinically reviewed and industry-listed as sharps safety device meeting O.S.H.A. requirements
- Packaged sterile
- Not made with natural rubber latex

Product No.		Description	Quantity	
Sterile	Non-Sterile		Sterile	Non-Sterile
D4506	D4506NS	#6 Safety Scalpel, Purple Slider	50/Cs	100/Cs
D4510	D4510NS	#10 Safety Scalpel, Black Slider	50/Cs	100/Cs
D4511	D4511NS	#11 Safety Scalpel, Blue Slider	50/Cs	100/Cs
D45E11	D45E11NS	#E11 Safety Scalpel, Orange Slider	50/Cs	100/Cs
D4514	D4514NS	#14 Safety Scalpel, Yellow Slider	50/Cs	100/Cs
D4515	D4515NS	#15 Safety Scalpel, Green Slider	50/Cs	100/Cs
D4515C	D4515CNS	#15C Safety Scalpel, Brown Slider	50/Cs	100/Cs
D4516	D4516NS	#16 Safety Scalpel, White Slider	50/Cs	100/Cs



* U.S. Patent # 5,531,754;
U.S. Patent # 6,254,621

BLADES & SCALPELS



** U.S. Patent # 5,531,754; U.S. Patent # 6,254,621*

#4 Safety Scalpels with Retractable Blades

- #4 handle houses larger stainless steel blades
- Patented* safety scalpel that allows a controlled retraction of the entire blade
- Permanent lock position for disposal on all versions
- Color-coded sliders for easy identification
- Packaged sterile
- Not made with natural rubber latex

Product No.		Description	Quantity	
Sterile	Non-Sterile		Sterile	Non-Sterile
D4520	D4520NS	#20 Safety Scalpel, Grey Slider	50/Cs	100/Cs
D4521	D4521NS	#21 Safety Scalpel, Lilac Slider	50/Cs	100/Cs
D4522	D4522NS	#22 Safety Scalpel, Sky Blue Slider	50/Cs	100/Cs
D4523	D4523NS	#23 Safety Scalpel, Lime Green Slider	50/Cs	100/Cs
D4524	D4524NS	#24 Safety Scalpel, Orange Slider	50/Cs	100/Cs

BLADES & SCALPELS

Swann-Morton® SCALPEL BLADES AND HANDLES



- Complete line of carbon steel and stainless surgical blades and quality stainless steel handles
- Manufactured in Sheffield, England, for DeRoyal by Swann-Morton®, the world leader in scalpel production
- Built on reputation for exceptionally high and consistent standards of durability and sharpness

Swann-Morton® Scalpel Handles

- Quality stainless steel
- Ideal weight and shape
- Slip-resistant grip
- Four sizes to accommodate a range of blade sizes
- Reusable
- Not made with natural rubber latex

Product No.	Description	Quantity
D6903	#3 Handle, Stainless Steel, Non-Sterile	40/Cs
D6904	#4 Handle, Stainless Steel, Non-Sterile	40/Cs
D0909	#9 Handle, Stainless Steel, Non-Sterile	40/Cs



BLADES & SCALPELS

Swann-Morton® FINE SCALPEL BLADES AND HANDLES



Fine Blades

- The Fine range has small, slim, precise blades in stainless steel for delicate operations
- Fine blades are used in disciplines such as reconstruction, ophthalmic, cardiology, and hair restoration
- Packaged sterile
- Not made with natural rubber latex

Product No.	Description	Quantity
D5901	Fine Surgical Blades, Stainless Steel, SM61	100/Cs
D5902	Fine Surgical Blades, Stainless Steel, SM62	100/Cs
D5903	Fine Surgical Blades, Stainless Steel, SM63	100/Cs
D5904	Fine Surgical Blades, Stainless Steel, SM64	100/Cs
D5905	Fine Surgical Blades, Stainless Steel, SM65	100/Cs
D5906	Fine Surgical Blades, Stainless Steel, SM64A	100/Cs
D5907	Fine Surgical Blades, Stainless Steel, SM67	100/Cs
D5908	Fine Surgical Blades, Stainless Steel, SM68	100/Cs
D5909	Fine Surgical Blades, Stainless Steel, SM69	100/Cs
D5911	Fine Surgical Blades, Stainless Steel, Singel Bevel, 61	100/Cs
D5912	Fine Surgical Blades, Stainless Steel, Singel Bevel, 62	100/Cs
D5921	Fine Surgical Blades, Stainless Steel, SP90	100/Cs
D5922	Fine Surgical Blades, Stainless Steel, SP91	100/Cs

BLADES & SCALPELS

Swann-Morton® SCALPEL BLADES

Swann-Morton® Scalpel Blades

- Available in carbon steel and stainless steel
- For use with #3 and #4 fitted handles
- All highest usage blade sizes available
- Individually foil wrapped
- Not made with natural rubber latex
- All packed 200/Case



Carbon Steel Blades

Product No.	Blade Size
D2201	SG3 Blade*
D6201	#10 Blade
D6203	#11 Blade
D6204	#12 Blade
D6218	#12D/12B Blade
D6205	#15 Blade
D6219	#15C Blade
D6206	#20 Blade
D6207	#21 Blade
D6208	#22 Blade
D6209	#22A Blade
D6210	#23 Blade
D6211	#24 Blade
D6215	#25A Blade

*Skin Graft Blade



Stainless Steel Blades

Product No.	Blade Size
D6301	#10 Blade
D6302	#10A Blade
D6303	#11 Blade
D6304	#12 Blade
D6318	#12D Blade
D6305	#15 Blade
D6321	#15C Blade
D6392	#15T Blade
D6306	#20 Blade
D6307	#21 Blade
D6308	#22 Blade
D6310	#23 Blade
D6311	#24 Blade

BLADES & SCALPELS

Swann-Morton® FINE SCALPEL BLADES AND HANDLES

SF13




SF23



SF3



SF4



Fine Handles

- Suited for encounters in such disciplines as reconstruction, ophthalmic, cardiology, and hair restoration
- Four special handles with independently screwed on top sections to facilitate blade-changing
- Packaged non-sterile
- Not made with natural rubber latex

Product No.	Description	Size
D6061	Fine Surgical Handle, SF13	10/Cs
D6062	Fine Surgical Handle, SF23	10/Cs
D6053	Fine Surgical Handle, SF3	10/Cs
D6054	Fine Surgical Handle, SF4	10/Cs


Disposable Scalpels

- Swann-Morton® scalpels with stainless steel blades
- Disposable affordability with high quality
- Packaged sterile
- Not made with natural rubber latex

Product No.	Description	Quantity
D6501	#10 Disposable Stainless Steel Scalpel	40/Cs
D6503	#11 Disposable Stainless Steel Scalpel	40/Cs
D2383	#11 Mini-Disposable Stainless Steel Scalpel	40/Cs
D6505	#15 Disposable Stainless Steel Scalpel	40/Cs
D6507	#21 Disposable Stainless Steel Scalpel	40/Cs
D6508	#22 Disposable Stainless Steel Scalpel	40/Cs
D6510	#23 Disposable Stainless Steel Scalpel	40/Cs
D6511	#24 Disposable Stainless Steel Scalpel	40/Cs



SURGICAL CONVENIENCE



Magnetic Instrument Pad

- Helps reduce accidental sharps injuries associated with instrument passing
- Available with or without white transfer zone for instrument placement
- Encourages hands-free passing of sharps as part of safety programs
- Not made with natural rubber latex

Without Transfer Zone

Product No.	Description	Quantity
25-001	Instrument Pad, Magnetic (10" x 16")	30 Pks/Cs
25-002	Instrument Pad, Magnetic (20" x 16")	20 Pks/Cs

With Transfer Zone

Product No.	Description	Quantity
25-003	Instrument Pad, Magnetic w/Transfer Zone (10" x 16")	30 Pks/Cs
25-004	Instrument Pad, Magnetic w/Transfer Zone (20" x 16")	20 Pks/Cs

Non-Magnetic Instrument Pad

- Provides sterile neutral field for the safe placement of instruments
- Instruments cling to special slip-proof surface
- Lies securely over drapes, yet light in weight
- May be used to cover mayo stand to reduce noise or rolled up on mayo for instrument positioning
- Not made with natural rubber latex

Product No.	Description	Quantity
25-011	Instrument Pad, 20" x 16"	30/Cs
25-012	Pink Instrument Pad, 12" x 16"	25/Cs
25-012NS	Non-Sterile Pink Instrument Pad, 12" x 16"	200/Cs



SURGICAL CONVENIENCE

Pre-Op Mini Marker

- Mini marker satisfies Joint Commission criteria for prevention of wrong site/wrong surgery
- Less waste than full size skin marker
- Reduces opportunity for cross contamination by lessening the temptation to reuse
- Comfortable thickness and length for ease of use
- Convenient storage – 100 each bin
- Non-tattooing ink, non-toxic ink
- Less expensive than standard skin marker
- Not made with natural rubber latex

Standard Ink	Description	Quantity
26-017NS	Pre-Op Skin Marker	800/Cs
26-018NS	Pre-Op Skin Marker	200/Cs
Prep-Resistant Ink	Description	Quantity
26-018PRNS	Pre-Op Skin Marker, Non-Sterile	200/Cs



Skin Markers and Rulers

- For smooth, clear line markings on patient's skin
- Prep-resistant gentian violet ink
- Non-toxic and non-tattooing
- Each barrel legibly ruled in centimeters
- Flexible 15cm ruler allows measurement of body curves
- Optional blank labels for identification of medications or specimens
- Individually packed in sterile peel pouches for fiber-free introduction into sterile field
- Dispenser boxes minimize storage and maximize convenience
- Not made with natural rubber latex

Standard Ink	Description	Quantity
26-001	Skin Marker, Standard Tip w/ Ruler (Blue Barrel)	50/Cs
26-002	Skin Marker, Fine Tip w/ Ruler (White Barrel)	50/Cs
26-003	Skin Marker, Standard Tip w/ Ruler & Labels (Blue Barrel)	50/Cs
26-004	Skin Marker, Fine Tip w/ Ruler & Labels (White Barrel)	50/Cs
26-005	Skin Marker, Standard Tip (Blue Barrel)	50/Cs
26-006	Skin Marker, Fine Tip (White Barrel)	50/Cs
32-345	Sterile 6" Ruler	30/Cs
Prep-Resistant Ink	Description	Quantity
26-001PR	Skin Marker, Standard Tip w/ Ruler	50/Cs
26-003PR	Skin Marker, Standard Tip w/ Ruler & Labels	50/Cs
26-005PR	Skin Marker, Standard Tip	50/Cs
26-005PRNS	Skin Marker, Standard Tip, Non-Sterile	100/Cs



Dual-Tip Skin Markers

- Ideal solution for writing medication labels and marking skin on the sterile field
- Supports Joint Commission (JCAHO) compliance for labeling medications/solutions on the sterile field
- Consolidates label marker and skin marker into one, making it economical and reducing waste
- Not made with natural rubber latex

Product No.	Description	Quantity
26-019	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip	50/Cs
26-019NS	Non-Sterile Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip	200/Cs
26-020	Dual Tip Skin & Utility Marker, Regular & Ultra Fine Tip w/Ruler & Labels	50/Cs



WARNING - This product can expose you to chemicals including gentian violet, which is known to the state of California to cause cancer. For more information, go to www.P65Warnings.ca.gov.


WOUND CLOSURE



Rotating Head Skin Stapler

- Rotating head offers complete 360° degree rotation for the perfect angle
- Ergonomic pistol-grip style design
- Utilizes patented "cat ear" staples with raised corners to:
 - Prevent rolling
 - Float above incision to reduce scarring cross marks
 - Make removal easier
- Not made with natural rubber latex
- Sterile

DISPOSABLE SKIN STAPLE INFORMATION FORMED STAPLE DIMENSIONS



CARTRIDGE CONFIGURATIONS
35 staples

Product No.	Description	Quantity
25-1001	Rotating Head 35-Count Skin Stapler	12/Cs, 6 Pks/Bx, 2 Bxs/Cs

WOUND CLOSURE



Premium Fixed Head Skin Stapler

- Fixed Head Stapler has clear viewing window for maximum visibility and 35-staple count
- All staplers have “pre-cocking” capacity for precision placement
- High quality, 316L stainless steel staples
- DeRoyal® Staples:
 - Maximize sharpness to reduce tissue trauma
 - Have wider tissue draw
 - Are taller for deeper tissue hold
- Not made with natural rubber latex
- Sterile



Product No.	Description	Quantity
25-3001*	Fixed Head 35-Count Skin Stapler, Wide Staple dimensions after closure: 4.7mm high x 6.5mm wide	24/Cs, 6 Pks/Bx, 4 Bxs/Cs
25-3001C	Premium Economy Fixed Head 35-Count Skin Stapler, Wide Staple dimensions after closure: 4.7mm high x 6.5mm wide	24/Cs, 6 Pks/Bx, 4 Bxs/Cs
25-3002*	Fixed Head 35-Count Skin Stapler, Regular Staple dimensions after closure: 3.7mm high x 5.7mm wide	24/Cs, 6 Pks/Bx, 4 Bxs/Cs
25-3002C	Premium Economy Fixed Head 35-Count Skin Stapler, Regular Staple dimensions after closure: 3.7mm high x 5.7mm wide	24/Cs, 6 Pks/Bx, 4 Bxs/Cs

⚠️ *WARNING - This product can expose you to chemicals including di(2-ethylhexyl)phthalate, which is known to the State of California to cause cancer and birth defects or other reproductive harm. For more information, go to www.P65Warnings.ca.gov.

WOUND CLOSURE

Deluxe Fixed Head Skin Stapler

- Contains a counter with no reversible function
- Open staple measures 11.8mm x 3.4mm
- Closed staple measures 6.8mm x 4.0mm
- Surgical-grade stainless steel staples
- Ergonomic handle design
- Open-ended for visibility and ease of placement and positioning
- Gear-down mechanism to slow staple formation and minimizes tearing during skin puncture
- Audible click to indicate staple is fully formed
- Soft and even trigger pull
- Easy-grasp packaging facilitates sterile presentation
- Not made with natural rubber latex
- Sterile



Product No.	Description	Quantity
25-3003	Deluxe 35-Count Skin Stapler	30/Cs

Skin Staple Extractor

- For the safe, easy removal of skin staples after wound healing
- Packaged sterile with all items needed for convenient aseptic staple removal, including alcohol prep pad and 2" x 2" gauze
- Not made with natural rubber latex
- Sterile

Product No.	Description	Quantity
TSR-1	Skin Staple Remover	12 Pks/Bx, 2 Bxs/Cs


* Apicare is a registered trademark of Apicare, Inc.



Curved upper jaw does not poke into patient.

Grooved jaws prevent staple from slipping during extraction.

SURGICAL ESSENTIALS



Sponge Counter Bag

- Allows continuous sponge counting and containment to minimize exposure to bloodborne pathogens
- Easy open pockets for the simple disposal of used sponges
- Accommodates five laparotomy sponges or ten 4" x 4"s
- Clear polyethylene pockets allow sponges to be viewed for visual estimate of patient blood loss
- Additional bags can be placed over filled ones as necessary
- Blue-back bags available with seal to contain sponges for safe disposal.
- Not made with natural rubber latex

Product No.	Description	Quantity
25-500	Sponge Counter Bags, Clear Back, Non-Sterile	50 Bags/Bx, 5 Bxs/Cs
25-551	Sponge Counter Bags, Blue Back w/Adhesive Seal, Non-Sterile	50 Bags/Bx, 5 Bxs/Cs
25-501	Mounting Bracket	1 each

SURGICAL ESSENTIALS

Instrument Pouch

- Keeps frequently used instruments out of the way yet in surgeon's reach as needed
- Comes with three pockets for multiple instruments
- Extra long sleeve securely holds complete laparoscopic instrument without dangling or danger of dropping
- Two surfaces of aggressive double-sided tape for secure placement
- Not made with natural rubber latex


Product No.	Description	Quantity
25-303	3-Pocket, 18.5" x 9.5"	30/Cs
25-302	2-Pocket, 11.5" x 7"	30/Cs



Cautery Tip Cleaner

- A convenient debris removal system
- Allows quick and easy cleaning of instruments
- Coated scratch pad eliminates disruptive scraping of electrodes
- Simple design applicable to mono or bipolar units
- Radiopaque element ensures x-ray detectability
- Not made with natural rubber latex

Product No.	Description	Quantity
32-200	2"x 2"	100 Pks/Cs
30-200	2"x 2"	240 Pks/Cs



Instrument Wipe

- For scratch-free cleaning of ophthalmic and other delicate instruments
- Non-linting
- Soft, sponge-like material (polyvinyl alcohol)
- Highly absorbent
- Not made with natural rubber latex

Product No.	Description	Quantity
25-201	Instrument Wipe, Polyvinyl Alcohol (2 1/2" x 2 1/2")	120/Cs



SAFETY DEVICES

Scalpel Holder

- Maintains scalpel in “blade-down” position to reduce the chance of injury
- Accommodates up to three scalpels
- Readily adheres to mayo stand or back table cover
- Holds any size scalpel handle and scalpel blade
- Not made with natural rubber latex

Product No.	Description	Quantity
30-920	Scalpel Holder	100 Pks/Cs



Safety Strainers

- Allows the removal of contaminated instruments from basins without requiring personnel to place hands in contaminated fluid
- Fits standard ring basins
- Does not interfere with normal ring basin use
- Not made with natural rubber latex

Product No.	Description	Quantity
89-5014	Rigid, Sterile Strainer	8/Cs
89-5020	Rigid, Strainer, Non-Sterile	100/Cs
30-910	Flexible, Sterile Strainer	50/Cs
30-910NS	Flexible Strainer, Non-Sterile	100/Cs



SAFETY DEVICES

Needle Cap Caddie™

- Device for needle recapping and exchange
- Accommodates all sizes and brands of hypodermic, spinal and I.V. needles
- Positions up to 10 needles at a time
- Indentations securely grip and station cap during needle removal and re-capping
- Helps to eliminate accidental needles sticks to clinician during recapping
- Adhesive backing for stable placement
- Patent pending
- Not made with natural rubber latex

Product No.	Description	Quantity
30-940	Sterile	30/Cs
30-940NS	Non-Sterile	30/Cs



HOW TO USE NEEDLE CAP CADDIE™



Step 1: Insert capped needle with syringe attached into appropriate hole in Needle Cap Caddie.



Step 2: Lift syringe with needle out of cap for use.



Step 3: Cap will stay in Needle Cap Caddie ready for safe, one-handed re-capping.



Step 4: To exchange needles, twist syringe to unscrew luer lock.



Step 5: Lift syringe away from cap containing the used needle and place another needle into Needle Cap Caddie.

HANDS-FREE PASSING



SharpStop™ Sharps Transfer Tray

- A disposable device to transfer and pass commonly used sharp instruments
- Use where controls are required to reduce employee exposure to bloodborne pathogens
- Helps to prevent sharps injuries during invasive procedures
- Facilitates compliance to OSHA Directive CPL 2-244D
- Accommodates wide range of scalpels, sutures, trocars, forceps, K-wires, syringes, and retractors
- Underside accommodates the practical positioning and safe passing of multiple instruments
- Designed by a scrub nurse and proven user-friendly for entire operating team
- Accommodates handle sizes 3,4,5,7, and 9 with loaded scalpel blades 10,11,12,15, 20, 22 and 23
- Center finger guide cradles scalpel for safe removal without looking away from surgical field
- Sloping sides automatically rotate and direct dropped scalpels into safety position
- Center guide accommodates up to size 9 gloved hand
- Specially designed central recess is wide enough to accommodate scalpel blades but does not allow fingers to touch blades
- Not made with natural rubber latex

Product No.	Description	Quantity
25-2244	Sharps Transfer Tray	24/Cs
25-2244NS	Sharps Transfer Tray, Non-Sterile	24/Cs
25-2244V	Sharps Transfer Tray w/Hook-and-Loop	24/Cs

SURGICAL CLIPS AND CLAMPS

Multi-Purpose Tube Holder



DeRoyal is pleased to offer a Multi-Purposed Tube Holder that frees clinicians from the task of securing lines with clumsy, make-shift tape constructions that are hard to secure and harder to keep in place and remove.

- Designed for quick and easy placement of all sizes of tubing and lines
- Easy-to-apply adhesive base attaches firmly to skin, cloth or drapes
- Gentle to skin and easily removed without painful pulling
- Holds lines securely in place once line is positioned on patient
- Secure placement for positioning to create slack without bladder traction, thus decreasing risk of urinary track infection
- Not made with natural rubber latex

For use with:


- Jejunostomy Tubing
- Leg Bag Tubing
- Subclavian Port Lines
- Wound drainage tubing
- Peritoneal dialysis exit site tubing
- A variety of other tubes and catheters

Product No.	Description	Quantity
30-1002	Tube Holder, Medium, Non-Sterile	50 Pks/Cs
30-1003	Tube Holder, Large, Non-Sterile	50 Pks/Cs




Competitors' Products

Other tube holders allow the tubing line to slide and dislodge from place.



DeRoyal® Tube Holder

DeRoyal® Tube Holder's adhesive surface holds tubing securely in place without slipping or displacement.



DeRoyal® Multi-Purpose Tube Holder can be attached to skin, drapes or cloth for use during a variety of procedures to accommodate numerous styles of lines.



SURGICAL CLIPS AND CLAMPS

Towel Clamp

- Positions drapes, towels and sponges in sterile field
- Built-in wire holder for cautery pencil or other wires
- Special bridge in curved jaws clamp and occlude a variety of tubes
- Strong and durable, the nylon reinforced plastic will not break
- Not made with natural rubber latex


Product No.	Description	Quantity
30-901	Towel Clamp	1/Pk, 80 Pks/Cs
30-902	Towel Clamp	2/Pk, 120 Pks/Cs
30-900NS	Towel Clamp, Non-Sterile	100/Cs
30-9010NS	Towel Clamp, Stainless	12/Cs
	Ball & Cup, 5 1/4", Non-Sterile	



Tube Positioner

- Allows stationary placement of tubes, lines, and cords in sterile field
- Extended adhesive peel-off backing for anchored placement
- Easily removed by gloved personnel
- Raypaque element ensures x-ray detectability
- Not made with natural rubber latex

Product No.	Description	Quantity
30-802	Tube/Cord Holder	25 Pks/Bx, 2 Bxs/Cs



SURGICAL CLIPS AND CLAMPS

Breathing Circuit Support Clamp

- Holds breathing circuits and suction tubes in place
- Relieves the “tug” on the ET (endotracheal tube) connection
- Reduces risk of kinks in ET tube, more comfortable for ICU patient
- Reduces chance of a disconnection
- Provides suction handle support
- Suction handle is more easily accessible to staff and patient
- Easy to use and holds tubing securely to patient's gown, bedding or drape
- Convenient support and easily positioned by staff
- Not made with natural rubber latex

Product No.	Description	Quantity
JLC-3NS	Breathing Circuit Support Clamp, Non-Sterile	100/Cs



Anesthesia Drape “Peach” Clip

- Used to attach surgical drape to IV poles to create barrier at patient's head
- Can be used for a variety of other clamping purposes
- Heavy duty holding power that will not slip
- Wide, easy-to-grip surface
- Non-rusting components
- Multiple colors
- Not made with natural rubber latex

Product No.	Description	Quantity
JPC-3NS	Assorted Colors, Non-Sterile	8/Bx








C-ARM DRAPES



C-Arm Drape

- Large, one-piece clear poly cover for complete coverage of C-Arm
- White poly adhesive straps secure drape in place
- Wide range of sizes available
- Packaged with telescoping fold
- Adhesive straps included to secure excess drape material
- Not made with natural rubber latex

Product #	Description	Quantity
30-0202	C-Arm Drape, 36" x 28" Rectangular Band Bag Style	20/Cs
		
30-0204	C-Arm Drape, 41" x 90" Mobile X-Ray Drape w/Extra Length	20/Cs
		
30-0201	C-Arm Drape, 41" x 74" Mobile X-Ray Drape	20/Cs
		

Product #	Description	Quantity
30-0209	C-Arm Drape, 41" x 140" Mobile X-Ray Drape	10/Cs
		
30-106M	C-Arm Drape, 41" x 84" Mini for OEC 6600	20/Cs
30-006M	C-Arm Drape, 41" x 84" Pack w/Keyboard & Footswitch Covers	20/Cs
		



Rubberband



Tape

CAMERA DRAPES



Endoscopic Camera & Laser Arm Drape

- For use in laparoscopic, arthroscopic and hysteroscopic procedures
- Provides sterile barrier
- Protects expensive equipment by eliminating need to sterilize camera
- Neatly fits around endoscope attached to camera
- Telescopic fold with cuff for easy draping and full cord coverage
- Allows visualization all the way through the drape's end
- Non-latex containing rubber bands and adhesive strips for secure closure
- Clear elasticized fenestration for secure fit on elasticized tip version
- Not made with natural rubber latex

Product


Description

Quantity

28-0401

Camera Drape, Laser Arm, 7" x 96"
Blue Elasticized Tip, 1/4" Opening


60/Cs



28-0402

Camera Drape, Laser Arm, 5" x 96"
Perforated Tip, 2" Opening

20/Cs



Tape



Elasticized



CASSETTE DRAPES

X-Ray Cassette Covers and Liners


- Clear poly bag with wide adhesive strip for closure
- Several sizes to fit a broad range of cassette widths and lengths
- Double Cassette Drape available with unique closure to ensure sterility
- Pull-away liner reduces chance of contamination once cassette is inserted
- Not made with natural rubber latex



Product #	Description	Quantity
30-0105	X-Ray Cassette Cover w/Adhesive Strip, 24" x 20", Standard	48/Cs
30-0101	X-Ray Cassette Cover w/Adhesive Strip, 24" x 25", Standard	20/Cs




30-0102	X-Ray Cassette Cover w/X-large Cuff, 21" x 40"	20/Cs
---------	--	-------




Step 1: Grasp both corners of drape to unfold. Position hands inside cuff at each corner in accepting position.




Step 2: Spread drape apart to allow accepted cassette to gently slide to bottom of drape.



Step 3: Fully extend drape.




Step 4: Fold drape until cuff fits opposite end of cassette.



Product # **Description**
 30-0103 X-Ray Cassette Cover w/Release Liner, 20" x 24"


Quantity

20/Cs



30-0104 X-Ray Cassette Cover w/Release Liner, 24" x 23"

20/Cs



Step 1: Position hands under white corner cuffs and spread drape apart to accept cassette.



Step 2: Insert cassette into drape through blue liner.



Step 3: Carefully remove liner by grasping outer sleeve and pulling up.



Step 4: Gently slide cassette to the bottom of the drape.



Step 5: Unfold white cuff.



Step 6: Fold cuff over cassette edge to close drape.

BAND BAGS

Fluoroscope/Banded Bags

- Clear polyethylene maintains sterility in Cath labs
- Ensures sterility and protects equipment during procedures
- A wide range of sizes available
- Can be used for many applications in Surgery, Special Procedure Labs, Labor & Delivery, G.I. Labs, Radiology, Emergency Department and others
- Rectangular, square and round configurations
- Not made with natural rubber latex



Product # Description

Rectangular Style




Product #	Description	Length	Width	Opening (Unstretched)	Quantity
77-103012T	Band Bag w/Rubberband & Tape	30"	12"	7"	25/Cs
77-103018T	Band Bag w/Rubberband & Tape	30"	18"	7"	25/Cs
77-103020T	Band Bag w/Rubberband & Tape	30"	20"	7"	25/Cs
77-103030T	Band Bag w/Rubberband & Tape	30"	30"	7"	25/Cs
77-103036T	Band Bag w/Rubberband & Tape	36"	30"	9"	25/Cs
77-103620T	Band Bag w/Rubberband & Tape	36"	20"	7"	25/Cs
77-103630T	Band Bag w/Rubberband & Tape	30"	36"	9"	25/Cs
77-103636T	Band Bag w/Rubberband & Tape	36"	36"	9"	25/Cs
77-103640T	Band Bag w/Rubberband & Tape	36"	48"	9"	25/Cs
77-104030T	Band Bag w/Rubberband & Tape	40"	30"	9"	25/Cs
77-104036T	Band Bag w/Rubberband & Tape	40"	36"	9"	25/Cs
77-104040T	Band Bag w/Rubberband & Tape	40"	40"	9"	25/Cs
77-104060T	Band Bag w/Rubberband & Tape	40"	60"	9"	25/Cs

Bouffant Style

Product #	Description	Depth	Quantity
77-100011	Band Bag, Round w/Sewn Elastic, Small	14"	25/Cs
77-100012	Band Bag, Round w/Sewn Elastic, Medium	22"	25/Cs
77-100013	Band Bag, Round w/Sewn Elastic, Large	27"	25/Cs
77-100025	Band Bag, Round w/Sewn Elastic, X-Large	50"	25/Cs

Fluoroscope Drape

77-100015T	Fluoroscope Drape Flat w/Tape	44"	35"	—	25/Cs
------------	-------------------------------	-----	-----	---	-------



Tape



Sewn Elastic

AUXILLARY DRAPES



Isolation Drape

- Protects intestines against trauma and contamination when outside the abdominal cavity during surgery
- Useful to isolate distal portions of extremities from contamination
- Other uses include the transport of specimens and to protect equipment during procedures and storage
- High quality, soft, flexible material that is strong and puncture-resistant
- Impervious vinyl film bag allows the visualization of organs at all times
- High tensile strength and elongation properties ensures the safety of the sterile environment without tearing
- Draw string to enclose organs to be kept moist throughout the duration of the case
- Moisture-proof surface prevents fluid loss through evaporation
- Not made with natural rubber latex

Product #	Description	Quantity
30-5510	Isolation Bag, Clear, 20" x 20"	40/Cs



Foot Switch Cover

- Non-sterile poly drawstring bag with two layers to protect O.R. foot controls from fluids
- Provides clarity for surgeon when operating foot controls
- Not made with natural rubber latex

Product #	Description	Quantity
30-0301NS	Foot Switch Cover, 16" x 18" x 3"	200/Cs

AUXILLARY DRAPES

Slush Drape

- Maintains sterile field
- Facilitates proper slush production and consistency
- Excellent and economical for cardiac and transplant procedures
- Not made with natural rubber latex

Product #	Description	Quantity
30-0602	Slush Drape, 54" x 54" (for use under sterile basin)	50/Cs



Case Cart Cover

- Non-sterile, 2 mil cover for critical items on case carts
- Clear poly perforated bag fits to size to keep items in place
- Not made with natural rubber latex

Product #	Description	Quantity
32-1385NS	Case Cart Cover, Clear, 48" x 72"	100/Roll
32-1386NS	Case Cart Cover, Opaque, 48" x 72"	100/Roll



Trademark Attributions

Aplicare is a registered trademark of Aplicare, Inc.

Swann-Morton and the Swann-Morton logo are registered trademarks of Swann-Morton Limited.

DeRoyal, the DeRoyal logo, Improving Care, Improving Business., DeBlade, Needle Cap Caddie, Saf-T-Glove, SharpStop, SurgiClick, and the Surgimate logo are trademarks or registered trademarks of DeRoyal Industries, Inc.



DeRoyal Industries, Inc.
200 DeBusk Lane, Powell, TN 37849 USA
800.DEROYAL or 865.938.7828
Customer Service: 800.251.9864
www.deroyal.com
Reprint #: 0-1433 | Rev. 7/19

©2019 DeRoyal Industries, Inc. All Rights Reserved.