

WE BRING EFFICIENCY TO YOUR PATIENT CARE

Fertility Treatment Range



We focus on patient safety by providing single-use procedure pack and contract manufacturing services that put you in control.



Clonallon

CUSTOMISING YOUR PATIENT CARE

Table of Contents

1. Who is Clonallon?
2. Why Choose a Customised Procedure Pack?
3. Embryo Transfer Pack
4. Oocyte Retrieval Pack
5. Oocyte Retrieval Pack Alternative
6. ET/IUI Pack
7. SIS/HyCoSy Pack
8. Procedure Pack
9. Premium Suture Pack
10. Gynae Exam Pack
11. Componentry and Instrumentation
12. Surgical Gowns





Clonallon

CUSTOMISING YOUR PATIENT CARE

Who is Clonallon?

Clonallon Laboratories Ltd was incorporated in 2001 to manufacture customised sterile surgical and medical procedure packs for global healthcare market.

Based in a modern facility in Northern Ireland, Clonallon manufacture sterile procedural packs in validated Class 7 Cleanrooms and all our processes and Quality Management Systems are accredited to ISO13485:2016 and MDD 93/42/EC. Distribution of Clonallon and third-party products is carried out from our warehouse and controlled by our own in-house bespoke software called the Management Information System. All Clonallon manufactured goods are provided with a Unique Device Indicator (UDI) which allows full traceability along the supply chain to an individual patient.

We focus on patient safety by providing single-use procedure packs and contract manufacturing services that put you in control. Firm clinical knowledge has allowed us to establish a trusted reputation in the healthcare industry and has enabled us to further develop and expand our services to client focused pack contract manufacturing. Our corporate ethos is defined by the intent to increase patient safety and reduce the overall burden of HCAs on our healthcare institutions.

Our products have been validated for a shelf-life of three years from the date of manufacture. We have worked closely with Fertility Clinics across the UK and Ireland to develop a product portfolio to suit their clinical requirements and you will find examples of these products within this brochure. However all of our procedure packs are fully customisable to suit your needs.

For any enquires please contact us for more information using the details below.



UK Office

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom

T: +44 (0)28 417 54531
E: michael@clonallon.com

Ireland Office

Blackthorn Business Park
Coe's Road, Dundalk
Co. Louth A91 EO3R
Ireland

T: +353 (0) 42 96 11202
E: michael@clonallon.com





Clonallon

CUSTOMISING YOUR PATIENT CARE

Why Choose a Customised Procedure Pack?

A Customised Procedure Pack delivers a specific set of disposable medical products provided in a single packed sterile tray format.

Benefits of using a Customised Procedure Pack:

1. Reduction of preparation times for procedures releasing expensive nursing time for other activities.
2. Standardise procedures and embed the prevention of Healthcare Acquired Infections in everyday practice.
3. Improve efficiency, reducing waiting times by increasing the number of procedures that can be completed
4. Reduced ordering and administration costs.
5. Release more time for direct patient care.
6. Caters for healthcare professional's preference as all packs are fully customisable.
7. Reduction of clinical waste in the theatre – environmentally friendly.
8. Complete traceability of all CE marked components used.
9. Eradication of errors in setting up for procedure.
10. Elimination of the risk of cross infection from a re-usable device.

Benefits of working with Clonallon:

1. Clonallon has been at the forefront of customised procedure pack design and manufacture for over 20 years. Our competency in this area is second to none.
2. We will work closely with you to develop your customised procedure pack to suit your clinical needs.
3. Samples are provided for evaluation.
4. Once the pack contents have been established, we will ensure optimum stock holdings are maintained to facilitate immediate delivery for recurring orders.
5. All our products are CE marked
6. Our Management Information System ensures complete traceability of components along the supply chain.
7. All our goods are fitted with a Unique Device Indicator.
8. All our packs are manufactured in our cleanroom facility in Northern Ireland. You can be confident of supply, quality and price competitiveness throughout the UK and Ireland.

*An example of the contents of one of our
Oocyte Retrieval Packs (Page 4)*





Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Embryo Transfer Pack

List of Components

- > 120ml Gallipot x 1
- > Transparent Kidney Dish, 800ml x 1
- > Rampley Sponge Holder 24cm x 1
- > Folded Hand Towel 450 x 490mm x 1
- > PE Wrap 60cm x 60cm x 1
- > Shallow Tray 25cm x 14cm x 2.5cm x 1
- > White Crepe 90cm x 90cm 60g x 1
- > White XRD Swabs 10cm x 10cm 32 ply Tied 5's x 5
- > Ball Non Woven Non Serile With Elastic Ring Size H x 5

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Oocyte Retrieval Pack

List of Components

- > 120ml Gallipot x 1
- > Transparent Kidney Dish, 800ml x 1
- > Rampley Sponge Holder 24cm x 1
- > Folded Hand Towel 450 x 490mm x 1
- > PE Wrap 60cm x 60cm x 1
- > Shallow Tray 25cm x 14cm x 2.5cm x 1
- > White Crepe 90cm x 90cm 60g x 1
- > White XRD Swabs 10cm x 10cm 32 ply Tied 5's x 5
- > Ball Non Woven Non Serile With Elastic Ring Size H x 5

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Oocyte Retrieval Pack Alternative

List of Components

- > Probe cover 42/125 x 610mm x 1
- > BIRR Needle Guide to fit GE E8C-RS Transducer x 1
- > White Crepe 90cm x 90cm 60g x 1

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

ET/UII Pack

List of Components

- > *Applicator Cotton Small 15cm x 2*
- > *White Crepe 90cm x 90cm 60g x 1*
- > *Green Flat Pouch 150x350mm x 1*
- > *Orchid Speculum Medium (Clear) x 1*

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

SIS/HyCoSy Pack

List of Components

- *Ramsey Sponge Holder 24cm x 1*
- *120ml Gallipot N/S x 1*
- *N/Woven White Swabs 7.5x7.5cm 4Ply x 4*
- *Tray Half Shallow 245x135x20mm x 1*
- *White Crepe 90cm x 90cm 60g x 1*
- *Non sterile Non woven Balls 40x45mm x 4*
- *Orchid Speculum Medium (Clear) x 1*

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- EO Sterilisation
- 3 Year Shelf Life
- Disposable / Single Use Only



Clonallon

CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Procedure Pack

List of Components

- > Rampley Sponge Holder 24cm x 1
- > 120ml Gallipot N/S x 1
- > White Crepe 90cm x 90cm 60g x 1

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Premium Suture Pack

List of Components

- > 120ml Gallipot x 2
- > Basin Emesis 700ml x 1
- > Utility Drape (Blue Tissue Laminate), 90cm x 90cm x 1
- > Folded Hand Towel 450 x 490mm x 1
- > White Gauze Plain Swab 10cm 32ply x 5
- > White Crepe 100cm x 100cm 60g x 1
- > Ball Non Woven Non Serile With Elastic Ring Size F x 5
- > Halstead Mosquito Forceps 12.5cm x 2
- > Forceps, Gillies 1x2 toothed, 15cm without pin x 1
- > Kilner Needle Holder 13cm x 1
- > S/S Iris Scissors straight x 1
- > Adson Forceps non tooth x 1
- > Iris Scissors cvd 11cm x 1
- > Shallow Tray 25cm x 14cm x 2.5cm x 1

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



Clonallon

CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202



www.clonallon.com



michael@clonallon.com



NSAI



Medical Devices
ISO 13485:2016
NSAI Certified



Quality Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Gynae Exam Pack

List of Components

- > Tray Half Deep 9.72x5.4x2ins x 1
- > Forceps Bonney Polypus and Ovum Forceps 24CM x 1
- > Rampley Sponge Holder 24cm x 1
- > Mayo Scissors cvd 17cm x 1
- > 60ml Gallipot graduated. N/S x 2
- > N/S Ball Non Woven Medium Diameter 4cm x 5
- > Tenaculum Forceps 1x1 Tooth x 1

Customisation

We can customise the contents of our procedure packs to suit your organisation's individual needs. We will work closely with you to provide optimal stock levels.

Quality Assurance

Clonallon Laboratories manufacture all procedure packs in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



 **Clonallon**
CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical
Devices
ISO 13485:2016
NSAI Certified



Quality
Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Componentry and Instrumentation

Individual Items Include:

- > *Ramsey Sponge Holder*
- > *Mayo Scissors*
- > *Iris Scissors*
- > *Halsted Mosquito Forceps*
- > *Forceps Vulsellum Teales Uterine Curved Toothed 2 : 3*
- > *Gillies Forceps*
- > *Cotton Applicators*
- > *Swabs*
- > *Kidney Dishes*
- > *Non Woven Gauze Balls*
- > *Ultrasound Machine Endocavity Probe Needle Guide*
- > *Bard Parker Handle Blades*
- > *Cusco Vaginal Speculum*
- > *Orchid Speculum*
- > *Kilner Needle Holder*
- > *Adson Toothed Forceps*
- > *Surgical Gowns (Non Sterile Also Available)*
- > *Ultrasonic Gel Sachets*

Please Contact Us For More Details On Our Complete Portfolio Of Available Components & Instruments



Quality Assurance

Clonallon Laboratories manufacture all medical devices in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > **EO Sterilisation**
- > **3 Year Shelf Life**
- > **Disposable / Single Use Only**



Clonallon

CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202



www.clonallon.com



michael@clonallon.com



NSAI



Medical Devices
ISO 13485:2016
NSAI Certified



Quality Management
ISO 9001:2015
NSAI Certified



Clonallon

CUSTOMISING YOUR PATIENT CARE

FERTILITY TREATMENT RANGE

Surgical Gowns



Colour Coded

To ensure the correct level of protection for your team all of our gowns and packaging are colour coded. This takes the guesswork out of gown selection.

Specifications

- + 40gsm SMS / 35gsm SMS (Lite)
- + Cotton Cuffs
- + Liquid Repellent
- + Breathable & Low Linting
- + Velcro Neck Fastening

Reinforced (Red Neckline)

High performance gowns for maximum protection

Standard (Green Neckline)

Standard Gowns Ideal for low fluid procedures

Standard Lite (Lite Green Neckline)

Basic Gowns Ideal for very low fluid procedures

Gown Range

- > Reinforced (Red Neckline)
- > Standard (Green Neckline)
- > Standard Lite (Lite Green Neckline)

Quality Assurance

Clonallon Laboratories manufacture all medical devices in strict accordance with ISO 13485 and the Medical Devices Directive 93/42/EEC.

Sterility

- > EO Sterilisation
- > 3 Year Shelf Life
- > Disposable / Single Use Only



Clonallon

CUSTOMISING YOUR PATIENT CARE

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom
Tel. +44 (0) 28 417 54531
Fax. +44 (0) 28 417 54563

Blackthorn Business Park
Coe's Road
Dundalk
A91 EO3R
Ireland
Tel. +353 (0) 42 96 11202

 www.clonallon.com  michael@clonallon.com



NSAI



Medical Devices
ISO 13485:2016
NSAI Certified



Quality Management
ISO 9001:2015
NSAI Certified

The Management Team



Dermot Dempster

Managing Director

Dermot is the Managing Director of Clonallon® Laboratories Ltd and holds a B.Sc. honours degree in Microbiology, Biochemistry and an MBA in International Business and has a special interest in developing and marketing medical devices which reduce or eliminate infections. Dermot has worked with Johnson & Johnson and with ICON Clinical Research a specialist clinical trial company. Dermot is a keen scuba diver in his spare time.



Michael Sands

Financial Director

Michael is the Financial Director of Clonallon® Laboratories Ltd and is responsible for the profitable running of Clonallon®. Michael is a member of the Chartered Accountants of Ireland (ACA) which was gained while working with ASM Accountants. Michael won Finance Director of the Year 2022 at the prestigious Irish Accountancy Awards. Michael is a keen golfer in his spare time.



Tiarnán Conway

Business Development Manager

Tiarnán holds a honours B.Sc. in Genetics and Genomics from Liverpool John Moores University and joined Clonallon® graduate entry programme. Tiarnán is responsible for the introduction of new products into existing and new markets. He also takes a leadership role in ensuring that tender submissions are submitted on a timely and correct basis.

Tiarnán is passionate about the role of new product development to fit the requirements of a modern surgical service. In his spare time Tiarnán is a keen sailor and instructor.

WE BRING EFFICIENCY TO YOUR PATIENT CARE

Fertility Treatment Range



UK Office

5 Milltown Industrial Estate
Warrenpoint
BT34 3FN
United Kingdom

T: +44 (0)28 417 54531
E: michael@clonallon.com

Ireland Office

Blackthorn Business Park
Coe's Road, Dundalk
Co. Louth A91 EO3R
Ireland

T: +353 (0) 42 96 11202
E: michael@clonallon.com

